


B-CARE ORAL CARE

Efficiency You Can See

MORE EFFECTIVE. EASY TO USE. COMPACT



MORE EFFECTIVE. EASY TO USE. COMPACT

B-care Oral Care is an innovative new approach to oral healthcare for intubated hospital patients. Featuring an electric toothbrush with integrated delivery of oral disinfectant or moisturizing solution, B-care Oral Care ensures a thorough clean and hydration of the whole oral cavity, so you can prevent infections and send your patients home with healthy teeth and gums.

THE NEED:

- Dental plaque biofilm can cause gingivitis and oral tissue bleeding.
- Oral biofilm is a significant risk factor for hospital acquired infections (HAI).
- During prolonged hospitalization, nursing staff are required to carry out a strict oral hygiene regime of brushing teeth and moisturizing inside the oral cavity.
- If oral care is not carried out effectively, the patient's general health may be impacted, even after they are released from hospital.

THE CHALLENGE:

- Manual toothbrushes that are currently used are less effective at cleaning teeth than an electric toothbrush.
- Common oral care protocols are difficult to implement and comply with.
- With common toothbrushes, it is difficult to reach and brush within the mouth of intubated patients.
- Most of the commonly provided disinfection fluids get wasted - oral disinfection and moisturizing fluids do not transfer well from external containers onto the brush and tend to drip down, meaning that less than the required amount is actually applied to the patient's mouth.

SONIC VIBRATING BRUSH



THE SOLUTION: B-CARE ORAL CARE

- Electric toothbrush vibrates 1000s of times per minute to clean teeth more effectively, with minimal hand movement.
- Unique cylindrical-shaped brush for ease of access, avoiding awkward hand twisting movements.
- Disinfectant or moisturizing fluid compartment contained within the brush handle compartment and dispensed through the brush tip, ensuring that it is all applied directly to the patient's mouth.
- Simple tracking of 24 hour oral hygiene protocol, with dials for counting number of brushing and moisturizing operations performed.
- Compact and flexible kit, adjustable to any hospital protocol and reducing required storage space.

